

Wedge Cushion Installation Guide

Have a question? Call 617-356-8106 or email support@gogrit.us.

Whether you'd like to increase/decrease dump angle, give your elbows some clearance from the wheels, give your legs some room, or if you just want some extra squish for your tush, the Wedge Cushion is an easy solution! We love these cushions because they address a number of comfort issues for a wide range of riders.

Option 1: Narrow edge forward

Position the narrow edge of the Wedge Cushion closest to the front of the seat pan. This configuration will lessen the dump angle and is often preferred by riders 5'4" or shorter.

Place the Wedge above or below the standard cushion, as desired.



Option 2: Wide edge forward

Position the wide edge of the Wedge Cushion closest to the front of the seat pan. This configuration will increase the dump angle and is often preferred by riders 6'3" or taller.

Place the Wedge above or below the standard cushion, as desired.



PRO TIP:

Is either cushion sliding? Place a few small strips of Velcro between it and the other cushion/seat pan —they'll keep everything in the right place.